

SPECIFICATIONS

2 large circular electrodes

• Connection: 2mm female banana plug

Diameter: 1.6 cmLength: 12 cmWeight: 25gNickel-free



Its inter-gluteal length prevents the risk of anorectal aspiration.

Its light weight means that patients are not inconvenienced during rehabilitation.

The wide circular electrodes in series are positioned to stimulate the entire anal musculature.

ration.
ed during rehabilitation.
timulate the entire

PRACTICAL BENEFITS OF THE PROBE

- 1 Electrical stimulation SEF
- Classical global stimulation of the anal sphincter muscles
- 2 Biofeedback BFB/EMG
- Recording of global electrical muscle activity in supine position
- Recording of global electrical muscle activity in orthostatic position

THERAPEUTIC APPLICATIONS OF THE PROBE

- 1 Electrical stimulation SEF
- Anal muscular insufficiency (direct)
 - → Phasic and/or tonic muscle strengthening
- Sphincter insufficiency
- → Bladder overactivity

2 - Biofeedback - BFB/EMG : lying or standing

- → Global anal awareness
 - Identification and measurement of anal muscular insufficiency
- Monitor progress of functional recovery
 - Control of reflex and voluntary locking during exercise
- 3 Treatment of incontinence due to anal muscular insufficiency
- Treatment of anal hypertonia