Technical sheet ANALYS BLUE

SPECIFICATIONS

· 2 circular electrodes

Connection: 2mm banana plug female or DIN 3 male

Diameter: 1.6 cmLength: 13.5 cm

Banana plug weight: 57gDIN 3 plug weight: 32g

PRESENTATION

Its inter-gluteal length prevents the risk of anorectal aspiration.

Its light weight means that patients are not inconvenienced during rehabilitation.

The circular electrodes are positioned over the anal sphincter.

The 3 "notches" make it easy to position and hold the probe over the anal sphincter.

PRACTICAL BENEFITS OF THE PROBE

- 1 Electrical stimulation SEF
- Classical R & L global stimulation of the anal sphincter muscles
- 2 Biofeedback BFB/EMG
- Recording of global electrical muscle activity in supine position
- Recording of global electrical muscle activity in orthostatic position

THERAPEUTIC APPLICATIONS OF THE PROBE

- 1 Electrical stimulation SEF
- Anal muscular insufficiency (direct)
 - → Phasic and/or tonic muscle strengthening
 - Sphincter deficiency

- → Bladder overactivity
- Trans-anal pubo-rectal reinforcement
- 2 Biofeedback BFB/EMG: lying or standing
 - → Global anal awareness
 - Identification and measurement of anal muscular insufficiency
- Monitor progress of functional recovery
 - Control of reflex and voluntary locking during exercise
- 3 Treatment of incontinence due to anal muscular insufficiency
- Treatment of anal hypertonia

